

Wachs-Umrechnungstabelle

Von Wachs- auf Legierungsgewicht. Wachsgewicht inkl. Gusskanalgewicht in Gramm ohne Berücksichtigung des Gusskegels.



Legierungen	g/cm³	0,2	0,3	0,4	0,5	0,6	0,7	0,8	0,9	1,0	1,2	1,4	1,6	1,8	2,0	2,2	2,4	2,6	2,8	3,0	3,2	3,4	3,6	3,8	4,0	4,2	4,4	4,6	4,8	5,0
Hochgoldhaltige Keramikaufbrennlegierungen																														
Bio KOOS-metic®	17,0	4	6	7	9	11	12	14	16	17	21	24	28	31	34	38	41	45	48	51	55	58	62	65	68	72	75	79	82	85
KOOS Biobest®	18,9	4	6	8	10	12	14	16	18	19	23	27	31	35	38	42	46	50	53	57	61	65	69	72	76	80	84	87	91	95
KOOS Bioextra®	19,1	4	6	8	10	12	14	16	18	20	23	27	31	35	39	43	46	50	54	58	62	65	69	73	77	81	85	88	92	96
KOOS Bio® 980	19,1	4	6	8	10	12	14	16	18	20	23	27	31	35	39	43	46	50	54	58	62	65	69	73	77	81	85	88	92	96
KOOS 960 plus	18,6	4	6	8	10	12	14	15	17	19	23	27	30	34	38	41	45	49	53	56	60	64	67	71	75	79	82	86	90	93
KOOS 950 extra	18,1	4	6	8	10	11	13	15	17	19	22	26	29	33	37	40	44	48	51	55	58	62	66	69	73	77	80	84	87	91
KOOS 950 supra	17,6	4	6	8	9	11	13	15	16	18	22	25	29	32	36	39	43	46	50	53	57	60	64	67	71	74	78	81	85	88
KOOS 900	17,8	4	6	8	9	11	13	15	17	18	22	25	29	33	36	40	43	47	50	54	57	61	65	68	72	75	79	82	86	89
Goldreduzierte Keramikaufbrennlegierungen																														
KOOS 530	13,9	3	5	6	7	9	10	12	13	14	17	20	23	26	28	31	34	37	39	42	45	48	51	53	56	59	62	64	67	70
KOOS 510	14,5	3	5	6	8	9	11	12	14	15	18	21	24	27	29	32	35	38	41	44	47	50	53	56	58	61	64	67	70	73
KOOS Eco-Classic	12,1	3	4	5	6	8	9	10	11	12	15	17	20	22	24	27	29	32	34	36	39	41	44	46	48	51	53	56	58	60
Palladiumbasis-Keramikaufbrennlegierungen																														
KOOS 400 M	12,0	3	4	5	6	8	9	10	11	12	15	17	20	22	24	27	29	32	34	36	39	41	44	46	48	51	53	56	58	60
KOOS 300	11,6	3	4	5	6	7	9	10	11	12	14	17	19	21	24	26	28	31	33	35	38	40	42	45	47	49	52	54	56	58
KOOS Unipal®	11,5	3	4	5	6	7	9	10	11	12	14	17	19	21	23	26	28	30	33	35	37	40	42	44	46	49	51	53	56	58
KOOS Eco-Palladent®	11,3	3	4	5	6	7	8	10	11	12	14	16	19	21	23	25	28	30	32	34	37	39	41	43	46	48	50	52	55	57
Universallegierungen (Guss- oder Aufbrennlegierungen)																														
KOOS Multiceralor® plus	16,7	4	6	7	9	11	12	14	16	17	21	24	27	31	34	37	41	44	47	51	54	57	61	64	67	71	74	77	81	84
KOOS Auronorm® plus	16,7	4	6	7	9	11	12	14	16	17	21	24	27	31	34	37	41	44	47	51	54	57	61	64	67	71	74	77	81	84
KOOS Austar® plus	16,7	4	6	7	9	11	12	14	16	17	21	24	27	31	34	37	41	44	47	51	54	57	61	64	67	71	74	77	81	84
KOOS Aurobest® supra	15,7	3	5	6	8	9	11	13	14	16	19	22	25	28	31	35	38	41	44	47	50	53	57	60	63	66	70	72	75	79
KOOS 69	15,6	4	5	7	8	10	11	13	15	16	19	22	25	29	32	35	38	41	44	47	50	54	57	60	63	66	69	72	75	78
KOOS 55	13,7	3	5	6	7	9	10	11	13	14	17	20	22	25	28	31	33	36	39	42	44	47	50	53	55	58	61	64	66	69
KOOS Ecobest®	13,9	3	4	6	7	8	10	11	13	14	17	19	22	25	28	31	33	36	39	42	44	47	50	53	56	58	61	64	67	70
KOOS Ecostar® M³	13,0	3	4	5	7	8	9	10	12	13	16	18	21	23	26	29	31	34	36	39	42	44	47	49	52	55	57	60	62	65
KOOS Ecodent®	10,6	3	4	5	6	7	8	9	10	11	13	15	17	20	22	24	26	28	30	32	34	37	39	41	43	45	47	49	51	53
Hochgoldhaltige Goldgusslegierungen																														
KOOS Classic Inlay plus	16,7	4	6	7	9	11	12	14	16	17	21	24	27	31	34	37	41	44	47	51	54	57	61	64	67	71	74	77	81	84
KOOS 72 plus	15,6	4	5	7	8	10	11	13	15	16	19	22	25	29	32	35	38	41	44	47	50	54	57	60	63	66	69	72	75	78
KOOS 71	15,3	4	5	7	8	10	11	13	14	16	19	22	25	28	31	34	37	40	43	46	49	53	56	59	62	65	68	71	74	77
KOOS 66	14,7	3	5	6	8	9	11	12	14	15	18	21	24	27	30	33	36	39	42	45	48	50	53	56	59	62	65	68	71	74
Goldreduzierte Goldgusslegierungen																														
KOOS 57	13,8	3	5	6	7	9	10	12	13	14	17	20	23	25	28	31	34	36	39	42	45	47	50	53	56	58	61	64	67	69